

CONOSTAN S-21 + K+Li+Sb 100 ppm

Art. ID CON-151-021-323
Unit 400 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,01	%		
	Aluminium (Al)	[7429-90-5]	100	µg/g		
	Antimony (Sb)	[7440-36-0]	0,01	%		
	Antimony (Sb)	[7440-36-0]	100	µg/g		
	Barium (Ba)	[7440-39-3]	0,01	%		
	Barium (Ba)	[7440-39-3]	100	µg/g		
	Boron (B)	[7440-42-8]	0,01	%		
	Boron (B)	[7440-42-8]	100	µg/g		
	Cadmium (Cd)	[7440-43-9]	0,01	%		
	Cadmium (Cd)	[7440-43-9]	100	µg/g		
	Calcium (Ca)	[7440-70-2]	0,01	%		
	Calcium (Ca)	[7440-70-2]	100	µg/g		
	Chromium (Cr)	[7440-47-3]	0,01	%		
	Chromium (Cr)	[7440-47-3]	100	µg/g		
	Copper (Cu)	[7440-50-8]	0,01	%		
	Copper (Cu)	[7440-50-8]	100	µg/g		
	Iron (Fe)	[7439-89-6]	0,01	%		
	Iron (Fe)	[7439-89-6]	100	µg/g		
	Lead (Pb)	[7439-92-1]	0,01	%		
	Lead (Pb)	[7439-92-1]	100	µg/g		
	Lithium (Li)	[7439-93-2]	0,01	%		
	Lithium (Li)	[7439-93-2]	100	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,01	%		
	Magnesium (Mg)	[7439-95-4]	100	µg/g		
	Manganese (Mn)	[7439-96-5]	0,01	%		
	Manganese (Mn)	[7439-96-5]	100	µg/g		
	Molybdenum (Mo)	[7439-98-7]	0,01	%		
	Molybdenum (Mo)	[7439-98-7]	100	µg/g		
	Nickel (Ni)	[7440-02-0]	0,01	%		
	Nickel (Ni)	[7440-02-0]	100	µg/g		
	Phosphorus (P)	[7723-14-0]	0,01	%		

Phosphorus (P)	[7723-14-0]	100	µg/g
Potassium (K)	[7440-09-7]	0,01	%
Potassium (K)	[7440-09-7]	100	µg/g
Silicon (Si)	[7440-21-3]	0,01	%
Silicon (Si)	[7440-21-3]	100	µg/g
Silver (Ag)	[7440-22-4]	0,01	%
Silver (Ag)	[7440-22-4]	100	µg/g
Sodium (Na)	[7440-23-5]	0,01	%
Sodium (Na)	[7440-23-5]	100	µg/g
Tin (Sn)	[7440-31-5]	0,01	%
Tin (Sn)	[7440-31-5]	100	µg/g
Titanium (Ti)	[7440-32-6]	0,01	%
Titanium (Ti)	[7440-32-6]	100	µg/g
Vanadium (V)	[7440-62-2]	0,01	%
Vanadium (V)	[7440-62-2]	100	µg/g
Zinc (Zn)	[7440-66-6]	0,01	%
Zinc (Zn)	[7440-66-6]	100	µg/g