

## **Custom Solution - PCB - Mix 11 - 10 µg/mL each of in n-Hexane**

Art. ID LM24-CUS-10252  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,4-Trichlorobenzene	[120-82-1]	10	µg/mL		
	Pentachlorobenzene	[608-93-5]	10	µg/mL		
	Hexachlorobenzene	[118-74-1]	10	µg/mL		
	Hexachloro-1,3-butadien e	[87-68-3]	10	µg/mL		
	2,4,4'-Trichlorobipheny l (PCB 28)	[7012-37-5]	10	µg/mL		
	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	10	µg/mL		
	2,2',4,5,5'-Pentachloro biphenyl (PCB 101)	[37680-73-2]	10	µg/mL		
	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	10	µg/mL		
	2,2',4,4',5,5'-Hexachlo robiphenyl (PCB 153)	[35065-27-1]	10	µg/mL		
	2,2',3,4,4',5'-Hexachlo robiphenyl (PCB 138)	[35065-28-2]	10	µg/mL		
	2,2',3,4,4',5,5'-Heptac hlorobiphenyl (PCB 180)	[35065-29-3]	10	µg/mL		