

Cast Iron, SUS, Ductile (Nodular) Iron Setting-Up Sample, Block 60 mm x 35 mm x 18 mm

Art. ID BAS-SUS-4/29
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Informational Value	Aluminium (Al)	[7429-90-5]	0,01	%		
Informational Value	Arsenic (As)	[7440-38-2]	<0,005	%		
Informational Value	Carbon (C)	[7440-44-0]	3,2	%		
Informational Value	Chromium (Cr)	[7440-47-3]	0,10	%		
Informational Value	Copper (Cu)	[7440-50-8]	0,78	%		
Informational Value	Magnesium (Mg)	[7439-95-4]	0,02	%		
Informational Value	Manganese (Mn)	[7439-96-5]	0,15	%		
Informational Value	Nickel (Ni)	[7440-02-0]	0,10	%		
Informational Value	Sulfur (S)	[7704-34-9]	0,01	%		
Informational Value	Antimony (Sb)	[7440-36-0]	<0,005	%		
Informational Value	Silicon (Si)	[7440-21-3]	2,7	%		
Informational Value	Tin (Sn)	[7440-31-5]	<0,005	%		
Informational Value	Titanium (Ti)	[7440-32-6]	0,06	%		
Informational Value	Vanadium (V)	[7440-62-2]	0,51	%		