

Custom Solution - 6 components in HNO3 5%

Art. ID LM24-CUS-71206-1x100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle), (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cadmium (Cd)	[7440-43-9]	100	µg/mL		
	Chromium (Cr)	[7440-47-3]	200	µg/mL		
	Copper (Cu)	[7440-50-8]	500	µg/mL		
	Manganese (Mn)	[7439-96-5]	500	µg/mL		
	Lead (Pb)	[7439-92-1]	200	µg/mL		
	Zinc (Zn)	[7440-66-6]	200	µg/mL		