

Domoic acid standard solution in Acetonitrile/water (1:19, v/v)

Art. ID	NRC-CRM-DA-h
Unit	0,5 mL
Deliverydetails	De Minimis UN1993 III E1

Description

NRC-CRM-DA-h is a certified calibration solution designed for analytical method development and accurate quantitation of DA. The concentration is suitable for preparing a dilution series for calibration of instruments such as liquid chromatography with detection by ultraviolet absorbance (LC–UV) or mass spectrometry (LC–MS), as well as for spiking control samples for recovery experiments.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Domoic acid (DA)	[14277-97-5]	97,1 ± 5,5	µg/mL		