

Custom Solution - 6 components in HNO3 2%

Art. ID LM24-CUS-71809-100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months. (4 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Sodium (Na) | [7440-23-5] | 1 | mg/L | | |
| | Tin (Sn) | [7440-31-5] | 5 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 0,02 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 0,01 | mg/L | | |
| | Aluminium (Al) | [7429-90-5] | 0,1 | mg/L | | |
| | Phosphorus (P) | [7723-14-0] | 100 | mg/L | | |