

Custom Solution - 6 components in HNO3 2%

Art. ID LM24-CUS-71809-100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months. (4 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na)	[7440-23-5]	1	mg/L		
	Tin (Sn)	[7440-31-5]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	0,02	mg/L		
	Cadmium (Cd)	[7440-43-9]	0,01	mg/L		
	Aluminium (Al)	[7429-90-5]	0,1	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		