

Cu-base solid

Art. ID DHJ-58A-Cu-Sn2E
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 97,8 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,014 | % | | |
| | Nickel (Ni) | [7440-02-0] | *0,04 | % | | |
| | Phosphorus (P) | [7723-14-0] | *0,07 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,0028 | % | | |
| | Silicon (Si) | [7440-21-3] | *0,002 | % | | |
| | Tin (Sn) | [7440-31-5] | 1,93 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,067 | % | | |
| | Bismuth (Bi) | [7440-69-9] | *0,0007 | % | | |
| | Antimony (Sb) | [7440-36-0] | *0,0004 | % | | |