

## **Cu-base solid**

Art. ID                      DHJ-58A-Cu-BMn40-1.5  
Unit                         disc  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Copper (Cu)       | [7440-50-8] | 55,9                | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,35                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 1,45                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | *40,9               | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | *0,0052             | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | *0,0015             | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,195               | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 0,0031              | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 0,0029              | %    |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 0,06                | %    |        |        |