

FAME Mix 37 components in n-Heptane

Art. ID	LM24-OP-4013014
Unit	1,5 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1206 II E2

Description

Validity: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methyl cis-10-heptadecanoate	[75190-82-8]	500	µg/mL		
	Methyl Butyrate	[623-42-7]	500	µg/mL		
	Hexanoic acid methyl ester	[106-70-7]	500	µg/mL		
	Octanoic acid-methyl ester	[111-11-5]	500	µg/mL		
	Decanoic acid methyl ester (Capric acid methyl ester)	[110-42-9]	500	µg/mL		
	Undecanoic acid methyl ester	[1731-86-8]	500	µg/mL		
	Dodecanoic acid methyl ester (Lauric acid methyl ester)	[111-82-0]	500	µg/mL		
	Methyl Tridecanoate	[1731-88-0]	500	µg/mL		
	Methyl myristoleate	[56219-06-8]	500	µg/mL		
	Tetradecanoic acid methyl ester (Myristic acid methyl ester)	[124-10-7]	500	µg/mL		
	cis-10-Pentadecenoic acid methyl ester	[90176-52-6]	500	µg/mL		
	Pentadecanoic acid methyl ester	[7132-64-1]	500	µg/mL		
	Methyl Palmitoleate	[1120-25-8]	500	µg/mL		
	Hexadecanoic acid methyl ester (Palmitic acid methyl ester)	[112-39-0]	1000	µg/mL		
	Heptadecanoic acid-methyl ester	[1731-92-6]	500	µg/mL		

yl ester				
Linolenic acid methyl ester	[301-00-8]	800		µg/mL
Linoleic acid-methyl ester	[112-63-0]	1000		µg/mL
Methyl Oleate	[112-62-9]	1000		µg/mL
Methyl linoelaidate	[2566-97-4]	500		µg/mL
9-trans-Octadecenoic acid methyl ester (Elaidic acid methyl ester)	[1937-62-8]	500		µg/mL
Methyl gamma-linolenate	[16326-32-2]	500		µg/mL
Octadecanoic acid methyl ester (Stearic acid methyl ester)	[112-61-8]	1000		µg/mL
5,8,11,14-Eicosatetraenoic acid methyl ester (Arachidonic acid methyl ester)	[2566-89-4]	800		µg/mL
cis-5,8,11,14,17-Eicosapentaenoic acid methyl ester	[2734-47-6]	800		µg/mL
cis-8,11,14-Eicosatrienoic acid methyl ester	[21061-10-9]	500		µg/mL
cis-11,14-Eicosadienoic acid methyl ester	[61012-46-2]	500		µg/mL
Methyl cis-11-eicosenoate	[2390-09-2]	500		µg/mL
Methyl eicosatrienoate	[55682-88-7]	500		µg/mL
Eicosanoic acid methyl ester (Arachidic acid methyl ester)	[1120-28-1]	500		µg/mL
Methyl Heneicosanoate	[6064-90-0]	500		µg/mL
Methyl-cis-4,7,10,13,16,19-Docosahexaenoate	[301-01-9]	800		µg/mL
cis-13,16-Docosadienoic acid methyl ester	[61012-47-3]	500		µg/mL
Methyl cis-13-docosenoate	[1120-34-9]	500		µg/mL

Docosanoic acid methyl ester (Behenic acid methyl ester)	[929-77-1]	500	µg/mL
Methyl Tricosanoate	[2433-97-8]	500	µg/mL
15-Tetracosanoic acid methyl ester (Nervonic acid methyl ester)	[2733-88-2]	500	µg/mL
Tetracosanoic acid methyl ester (Lignoceric acid methyl ester)	[2442-49-1]	500	µg/mL