

Okadaic acid, approx. 12 µg/mL, (15.6 ± 0.9) µg/g in Methanol (ISO 17034), Lot: 22-001

Art. ID CFG-CRM-02-OA
Unit approx. 0,5 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------|---------------|---------------------|------|--------|--------|
| | Okadaic acid sodium sal t | [209266-80-8] | | | | |