

Low Phosphorus Iron, SUS, Block 60 mm x 35 mm x 18 mm

Art. ID BAS-SUS-1/20
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,2	%		
	Silicon (Si)	[7440-21-3]	2,7	%		
	Manganese (Mn)	[7439-96-5]	0,45	%		
	Phosphorus (P)	[7723-14-0]	0,06	%		
	Sulfur (S)	[7704-34-9]	0,06	%		
	Chromium (Cr)	[7440-47-3]	0,48	%		
	Molybdenum (Mo)	[7439-98-7]	0,34	%		
	Nickel (Ni)	[7440-02-0]	0,21	%		
	Aluminium (Al)	[7429-90-5]	0,02	%		
	Copper (Cu)	[7440-50-8]	0,53	%		
	Tin (Sn)	[7440-31-5]	0,05	%		
	Titanium (Ti)	[7440-32-6]	0,013	%		
	Vanadium (V)	[7440-62-2]	0,04	%		