

Austentic (Ni-Resist), Cast Iron, RM, Block 48 mm x 42 mm x 12 mm

Art. ID BAS-NIRM 1/1
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	1,83 ± 0,02	%		
	Silicon (Si)	[7440-21-3]	3,26 ± 0,01	%		
	Manganese (Mn)	[7439-96-5]	6,74 ± 0,04	%		
	Phosphorus (P)	[7723-14-0]	0,058 ± 0,001	%		
	Silicon (Si)	[7440-21-3]	0,015 ± 0,001	%		
	Chromium (Cr)	[7440-47-3]	0,3 ± 0,005	%		
	Nickel (Ni)	[7440-02-0]	11,8 ± 0,2	%		
	Copper (Cu)	[7440-50-8]	0,2 ± 0,01	%		
	Cerium (Ce)	[7440-45-1]	0,022 ± 0,002	%		
	Magnesium (Mg)	[7439-95-4]	0,021 ± 0,001	%		