

Custom Solution - 3 components 100 mg/L each of in HNO3 2%

Art. ID LM24-CUS-72005
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months. (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		