

## **Siloxanes - Mix 8, 100 µg/mL in Acetone**

Art. ID LM24-N-13570-100AC5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hexamethylcyclotrisiloxane	[541-05-9]	100	µg/mL		
	Octamethylcyclotetrasiloxane	[556-67-2]	100	µg/mL		
	Decamethylcyclopentasiloxane	[541-02-6]	100	µg/mL		
	Dodecamethylcyclohexasiloxane	[540-97-6]	100	µg/mL		
	Hexamethyldisiloxane	[107-46-0]	100	µg/mL		
	Octamethyltrisiloxane	[107-51-7]	100	µg/mL		
	Decamethyltetrasiloxane	[141-62-8]	100	µg/mL		
	Dodecamethylpentasiloxane	[141-63-9]	100	µg/mL		