

## Organotin - Mix 8 (DIN EN ISO 17353:2004), 100 µg/mL in Methanol (based on active substance)

Art. ID LM24-N-13970-W100ME10  
 Unit 10 mL  
 Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substance	n-Butyltin trichloride	[1118-46-3]	100	µg/mL		
based on active substance	Dibutyltin-dichloride	[683-18-1]	100	µg/mL		
based on active substance	Tributyltin chloride	[1461-22-9]	100	µg/mL		
based on active substance	Tetrabutyltin	[1461-25-2]	100	µg/mL		
based on active substance	n-Octyltin-trichloride	[3091-25-6]	100	µg/mL		
based on active substance	Di-n-octyltin-dichloride	[3542-36-7]	100	µg/mL		
based on active substance	Triphenyltin chloride	[639-58-7]	100	µg/mL		
based on active substance	Tricyclohexyltin chloride	[3091-32-5]	100	µg/mL		