

Siloxanes - Mix 8, 50 µg/mL in Acetone

Art. ID LM24-N-13570-50AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hexamethylcyclotrisiloxane	[541-05-9]	50	µg/mL		
	Octamethylcyclotetrasiloxane	[556-67-2]	50	µg/mL		
	Decamethylcyclopentasiloxane	[541-02-6]	50	µg/mL		
	Dodecamethylcyclohexasiloxane	[540-97-6]	50	µg/mL		
	Hexamethyldisiloxane	[107-46-0]	50	µg/mL		
	Octamethyltrisiloxane	[107-51-7]	50	µg/mL		
	Decamethyltetrasiloxane	[141-62-8]	50	µg/mL		
	Dodecamethylpentasiloxane	[141-63-9]	50	µg/mL		