

Organotin - Mix 4 (DIN EN ISO 17353:2004), 1000 µg/mL in Methanol (based on active substance)

Art. ID LM24-N-13980-W1000ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substance	Mono-n-heptyltin trichloride	[59344-47-7]	1000	µg/mL		
based on active substance	Di-n-heptyltin-dichloride	[74340-12-8]	1000	µg/mL		
based on active substance	Tri-n-propyltin-chloride	[2279-76-7]	1000	µg/mL		
based on active substance	Tetra-n-propyltin	[2176-98-9]	1000	µg/mL		