

Heterocycles - Mix 20, 50 µg/mL in Methanol

Art. ID LM24-N-14911-50ME5

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acridone	[578-95-0]	50	µg/mL		
	Benzo(b)furan	[271-89-6]	50	µg/mL		
	1-Benzothiophen	[95-15-8]	50	µg/mL		
	Carbazole	[86-74-8]	50	µg/mL		
	Quinoline	[91-22-5]	50	µg/mL		
	2-Hydroxyquinoline	[59-31-4]	50	µg/mL		
	Dibenzofuran	[132-64-9]	50	µg/mL		
	2,3-Dimethylbenzofuran	[3782-00-1]	50	µg/mL		
	2,4-Dimethylquinoline	[1198-37-4]	50	µg/mL		
	2-Phenylphenol	[90-43-7]	50	µg/mL		
	Isoquinoline	[119-65-3]	50	µg/mL		
	1-Hydroxyisoquinoline	[491-30-5]	50	µg/mL		
	2-Methylbenzofuran	[4265-25-2]	50	µg/mL		
	3-Methylbenzothiophene	[1455-18-1]	50	µg/mL		
	2-Methylquinoline (Quin aldine)	[91-63-4]	50	µg/mL		
	2-Hydroxy-4-methylquino line	[607-66-9]	50	µg/mL		
	2-Methyldibenzofuran	[7320-51-6]	50	µg/mL		
	1-Methylisoquinoline	[1721-93-3]	50	µg/mL		
	Phenanthridinone	[1015-89-0]	50	µg/mL		
	Xanthone	[90-47-1]	50	µg/mL		