

## **Thiacloprid-d4 100 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-5690D-100AN1  
Unit 1 mL  
Deliverydetails De Minimis UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number     | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|----------------|---------------------|-------|--------|--------|
|                  | Thiacloprid-d4    | [1793071-39-2] | 100                 | µg/mL |        |        |