

Gabapentin-d4 10 µg/mL in Acetonitrile

Art. ID LM24-N-18000-0496D-10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|----------------|---------------------|-------|--------|--------|
| | Gabapentin-d4 | [1185039-20-6] | 10 | µg/mL | | |