

## **Ethyl caproate-d11 100 µg/mL in Methanol**

Art. ID LM24-N-14600-0226D-100ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ethyl Caproate-d11	[2159-19-5]	100	µg/mL		