

Sweeteners - Mix 6, 100 µg/mL in Methanol (based on active substance)

Art. ID LM24-N-R02191-W100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substance	Acesulfame potassium	[55589-62-3]	100	µg/mL		
ce						
based on active substance	Saccharin	[81-07-2]	100	µg/mL		
ce						
based on active substance	Sucralose	[56038-13-2]	100	µg/mL		
ce						
based on active substance	1H-Benzotriazole	[95-14-7]	100	µg/mL		
ce						
based on active substance	4-Methyl-1H-benzotriazole	[29878-31-7]	100	µg/mL		
ce	le					
based on active substance	5-Methyl-1H-benzotriazole	[136-85-6]	100	µg/mL		
ce	le					