

Tetradecamethylcycloheptasiloxane 100 µg/mL in Acetone

Art. ID LM24-N-13560-0085-100AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tetradecamethylcyclohep tasiloxane	[107-50-6]	100	µg/mL		