

## **Fenpropimorph-d3 10 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-5585D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fenpropimorph-d3	[1292815-71-4]	10	µg/mL		