

### **1,3,5-Trimethylbenzene-d12 100 µg/mL in Methanol**

Art. ID LM24-N-13400-3901D-100ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mesitylene-d12	[69441-16-3]	100	µg/mL		