

### **Formamide-d3 1000 µg/mL in Methanol**

Art. ID LM24-N-13700-0230D-1000ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Formamide-d3	[43380-64-9]	1000	µg/mL		