

## **2-Chloroethylether-d8 100 µg/mL in Methanol**

Art. ID LM24-N-10900-0166D-100ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bis(2-chloroethyl)ether -d8	[93952-02-4]	100	µg/mL		