

### **Bis-(2-chloroethoxy)-methane-d8 100 µg/mL in Methanol**

Art. ID LM24-N-13600-0938-100ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	bis(2-Chloroethoxy)-d8- methane	[93966-78-0]	100	µg/mL		