

4-Acetamidoantipyrine-d3 10 µg/mL in Acetonitrile

Art. ID LM24-N-18000-2295D-10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------|---------------|---------------------|-------|--------|--------|
| | 4-Acetamido Antipyrine- d3 | [342821-66-3] | 10 | µg/mL | | |