

**Bisoprolol-d7 hemifumarate 10 µg/mL in Methanol (based on active substance)**

Art. ID LM24-N-18000-2304D-W10ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substance	Bisoprolol-d7 Hemifumarate	[1310012-16-8]	10	µg/mL		