

Bisoprolol-d7 hemifumarate 100 µg/mL in Methanol (based on active substance)

Art. ID LM24-N-18000-2304D-W100ME1
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substance	Bisoprolol-d7 Hemifumarate	[1310012-16-8]	100	µg/mL		