

**(S)-(-)-Valsartan-d8 10 µg/mL in Acetonitrile**

Art. ID LM24-N-18000-1256D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

**Description**

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(S)-(-)-Valsartan-D8 (V aline-D8)	[1089736-72-0]	10	µg/mL		