

1,1,1-Trichloroethane-d3 100 µg/mL in Methanol

Art. ID LM24-N-10900-0771D-100ME2
Unit 2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,1,1-Trichloroethane-2 ,2,2-d3	[2747-58-2]	100	µg/mL		