

## **2,4'-DDT-d8 10 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-1301D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4'-DDT D8	[221899-88-3]	10	µg/mL		