

2,4-Dichlorotoluene-3,5,6-d3 100 µg/mL in Methanol

Art. ID LM24-N-13400-4479D-100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4-Dichlorotoluene-3,5 ,6-D3	[1398065-47-8]	100	µg/mL		