

2-Ethanolamine-d4 100 µg/mL in Methanol

Art. ID LM24-N-16100-0181D-100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ethanol-1,1,2,2-d4-amin e	[85047-08-1]	100	µg/mL		