

### **4-Acetamidoantipyrine-d3 10 µg/mL in Methanol**

Art. ID LM24-N-18000-2295D-10ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Acetamido Antipyrine- d3	[342821-66-3]	10	µg/mL		