

5a-Androstenone-d4 10 µg/mL in Acetonitrile

Art. ID LM24-N-18000-0059D-10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|------------|---------------------|-------|--------|--------|
| | 5a-Androstenone-d4 | [n/a] | 10 | µg/mL | | |