

Bisoprolol-d7 hemifumarate 10 µg/mL in Methanol

Art. ID LM24-N-18000-2304D-10ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisoprolol-d7 Hemifumarate	[1310012-16-8]	10	µg/mL		