

Clofibric acid-d4 10 µg/mL in Methanol

Art. ID LM24-N-18000-0273D-10ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Clofibric-d4 Acid	[1184991-14-7]	10	µg/mL		