

### **DL-Phenylalanine-d8 10 µg/mL in Methanol**

Art. ID LM24-N-18000-1656D-10ME10  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-Phenylalanine-d8	[29909-00-0]	10	µg/mL		