

## **DL-Phenylalanine-d8 100 µg/mL in Acetonitrile**

Art. ID LM24-N-18000-1656D-100AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-Phenylalanine-d8	[29909-00-0]	100	µg/mL		