

DNOC-d3 10 µg/mL in Cyclohexane

Art. ID LM24-N-15900-1921D-10CY10
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1145 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,6-Dinitro-2-Methyl-D3 -Phenol	[1219804-69-9]	10	µg/mL		