

DNOC-d3 10 µg/mL in Cyclohexane

Art. ID LM24-N-15900-1921D-10CY10
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1145 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------------|----------------|---------------------|-------|--------|--------|
| | 4,6-Dinitro-2-Methyl-D3 -Phenol | [1219804-69-9] | 10 | µg/mL | | |