

## **DNOC-d3 100 µg/mL in Acetone**

Art. ID LM24-N-15900-1921D-100AC5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,6-Dinitro-2-Methyl-D3 -Phenol	[1219804-69-9]	100	µg/mL		