

Isoprothiolane-d4 10 µg/mL in Acetone

Art. ID LM24-N-15900-2926D-10AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|----------------|---------------------|-------|--------|--------|
| | Isoprothiolane-d4 | [1715020-82-8] | 10 | µg/mL | | |