

N-Nitroso-N-methyl-4-aminobutyric acid-d3 100 µg/mL in Methanol

Art. ID LM24-N-16000-0069D-100ME1
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Nitroso-N-methyl-4-am inobutyric Acid-d3	[1184996-41-5]	100	µg/mL		