

Trichloroethylene-d1 100 µg/mL in Methanol

Art. ID LM24-N-10900-0801D-100ME2
Unit 2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Trichloroethylene-d	[13291-68-4]	100	µg/mL		