

### **Valsartan-d3 10 µg/mL in Acetonitrile**

Art. ID LM24-N-18000-2363D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

#### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Valsartan-d3	[1331908-02-1]	10	µg/mL		